

Anti-Rabggta Polyclonal Antibody For Wb, If/Icc, Elisa - Q92696

Item Number: CM0024484

Introduction

Rabbit polyclonal antibody against human RABGGTA (Rab geranylgeranyltransferase subunit alpha) validated for Western blot, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for prenylation studies.

[Learn More](#)

Parameter	Value
Product Name	RABGGTA Rabbit pAb
Remarks/Alias	PTAR3; RABGGTA
Species	Human
Gene ID (Human)	5875
Gene ID	5875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-360 of human RABGGTA (NP_878256.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92696
Protein Weight	65kDa
Shipping	Ice bag