

Anti-Rabggtb Polyclonal Antibody For Wb, Elisa - P53611

Item Number: CM0019057



Introduction

Rabbit polyclonal antibody against human RABGGTB (geranylgeranyl transferase type-2 subunit beta). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for prenylation and intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	RABGGTB Rabbit pAb
Remarks/Alias	GGTB; RABGGTB
Species	Human
GeneID (Human)	5876
GeneID	5876
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-331 of human RABGGTB (NP_004573.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53611
Protein Weight	37kDa
Shipping	Ice bag