

Anti-Rab13 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Hyi8

Item Number: CM0018719

Introduction

Rabbit polyclonal antibody against human RABL3 (PNCA5), validated for WB, IF/ICC, and ELISA. Reacts with human and mouse samples. Targets a small GTPase involved in KRAS signaling and ciliogenesis.

[Learn More](#)

Parameter	Value
Product Name	RABL3 Rabbit pAb
Remarks/Alias	PNCA5; RABL3
Species	Human
Gene ID	285282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-236 of human RABL3 (NP_776186.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5HYI8
Protein Weight	26kDa
Shipping	Ice bag