


Anti-Rac Gtpase-Activating Protein 1 (Rgap1) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H0H5

Item Number: CM0027125



Introduction

Rabbit polyclonal antibody against human RGAP1 (Rac GTPase-activating protein 1), validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, and rat RGAP1; predicted molecular weight 71 kDa. Suitable for cell cycle and cytokinesis research.

[Learn More](#)

Parameter	Value
Product Name	RGAP1 Rabbit pAb
Remarks/Alias	CYK4; CDAN3B; ID-GAP; HsCYK-4; MgcRacGAP; RGAP1
Species	Human
Gene ID (Human)	29127
Gene ID	29127
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 353-632 of human RGAP1 (NP_001119576.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0H5
Protein Weight	71kDa
Shipping	Ice bag