


Anti-Rac1/2/3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P63000 / P15153 / P60763

Item Number: CM0010315



Introduction

High-affinity rabbit monoclonal antibody against Rac1/2/3, suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat Rac proteins. Recombinant immunogen based on human RAC1 N-terminal region (P63000). Ideal for studies on GTPase signaling and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Rac1/2/3 Rabbit mAb
Remarks/Alias	MIG5; MRD48; Rac-1; TC-25; p21-Rac1; Rac1/2/3
Species	Human
GeneID (Human)	5879
GeneID	5879/5880/5881
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human RAC1 (P63000).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63000/P15153/P60763
Protein Weight	21kDa
Shipping	Ice bag