


Anti-Rac1/2/3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P63000 / P15153 / P60763

Item Number: CM0028674



Introduction

Anti-Rac1/2/3 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets Rac1, Rac2, Rac3 (P63000/P15153/P60763). Ideal for cell migration and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Rac1/2/3 Rabbit pAb
Remarks/Alias	MIG5; MRD48; Rac-1; TC-25; p21-Rac1; Rac1/2/3
Species	Human
GeneID (Human)	5879/ 5880/ 5881
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human Rac1/2/3 (NP_008839.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63000/P15153/P60763
Protein Weight	21kDa
Shipping	Ice bag