

Anti-Rac2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15153

Item Number: CM0015565



Introduction

Rabbit polyclonal antibody targeting human RAC2 (p21-Rac2), for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Detects endogenous RAC2 in hematopoietic cells; ideal for NADPH oxidase, phagocytosis, and GTPase signaling research.

[Learn More](#)

Parameter	Value
Product Name	RAC2 Rabbit pAb
Remarks/Alias	Gx; EN-7; IMD73A; IMD73B; IMD73C; HSPC022; p21-Rac2; RAC2
Species	Human
GeneID (Human)	5880
GeneID	5880
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-192 of human RAC2 (NP_002863.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15153
Protein Weight	21kDa
Shipping	Ice bag