

Anti-Rack1 Monoclonal Antibody For Wb And Elisa - P63244

Item Number: CM0010972



Introduction

Rabbit monoclonal antibody against human RACK1 (P63244), validated for WB and ELISA. Reacts with human, mouse, rat. Useful for translation, signaling, and cell cycle research.

[Learn More](#)

Parameter	Specification
Product Name	RACK1 Rabbit mAb
Remarks/Alias	H12.3; HLC-7; PIG21; GNB2L1; Gnb2-rs1; RACK1
Species	Human
Gene ID (Human)	10399
Gene ID	10399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-310 of human RACK1 (P63244).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63244
Protein Weight	35kDa
Shipping	Ice bag