

Anti-Rack1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P63244

Item Number: CM0019032



Introduction

Rabbit polyclonal antibody targeting human RACK1 (GNB2L1), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for investigating signaling, translation, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	RACK1 Rabbit pAb
Remarks/Alias	H12.3; HLC-7; PIG21; GNB2L1; Gnb2-rs1; RACK1
Species	Human
GeneID (Human)	10399
GeneID	10399
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-290 of human RACK1 (NP_006089.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63244
Protein Weight	35 kDa
Shipping	Ice bag