


Anti-Rad18 Monoclonal Antibody For Wb And Elisa - Q9Ns91

Item Number: CM0004112



Introduction

High-quality rabbit monoclonal antibody targeting Rad18 for WB and ELISA applications. Cross-reacts with human, mouse, and rat. Suitable for postreplication repair and DNA damage studies.

[Learn More](#)

Parameter	Value
Product Name	Rad18 Rabbit mAb
Remarks/Alias	RNF73; Rad18
Species	Human
GeneID (Human)	56852
GeneID	56852
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Rad18 (Q9NS91)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NS91
Protein Weight	56 kDa
Shipping	Ice bag