


Anti-Rad18 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ns91

Item Number: CM0027121



Introduction

CamelBio offers an Anti-RAD18 rabbit polyclonal antibody validated for Western blot and ELISA, targeting human RAD18 (RNF73) E3 ubiquitin-protein ligase. Essential for postreplication repair and DNA damage research studies.

[Learn More](#)

Parameter	Description
Product Name	Rad18 Rabbit pAb
Remarks/Alias	RNF73; Rad18
Species	Human
Gene ID (Human)	56852
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-495 of human Rad18 (NP_064550.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NS91
Protein Weight	56kDa
Shipping	Ice bag