

Anti-Rad21 Rabbit Monoclonal Antibody For Wb, Elisa, Chip, Chip-Seq - O60216

Item Number: CM0002297



Introduction

Recombinant rabbit monoclonal antibody against human Rad21, validated for WB, ELISA, ChIP, and ChIP-seq. Cross-reacts with mouse and rat. Ideal for cohesin, cell cycle, and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Rad21 Rabbit mAb
Remarks/Alias	MGS; HR21; MCD1; NXP1; SCC1; CDLS4; hHR21; HRAD21; Rad21
Species	Human
GeneID	5885
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 532-631 of human Rad21 (O60216).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat
SWISS	O60216
Protein Weight	72kDa
Shipping	Ice bag