

Anti-Rad21 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - O60216

Item Number: CM0022492



Introduction

Rabbit polyclonal antibody targeting human Rad21 (SCC1 cohesin subunit). Validated for WB, IHC-P, IF/ICC, IP, ELISA, ChIP. Cross-reacts with mouse and rat. Suitable for chromatin biology, cell cycle, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Rad21 Rabbit pAb
Remarks/Alias	MGS; HR21; MCD1; NXP1; SCC1; CDLS4; hHR21; HRAD21; Rad21
Species	Human
GeneID (Human)	5885
GeneID	5885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 281-420 of human Rad21 (NP_006256.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	O60216
Protein Weight	72kDa
Shipping	Ice bag