

Anti-Rad23A Monoclonal Antibody For Wb, If/Icc, Elisa - P54725

Item Number: CM0008889



Introduction

Rad23A Rabbit mAb validated for WB, IF/ICC, ELISA. Recognizes human and mouse Rad23A (HR23A). Research-grade IVD raw material for DNA repair, histone demethylation, and proteasomal degradation studies.

[Learn More](#)

Parameter	Value
Product Name	Rad23A Rabbit mAb
Remarks/Alias	HR23A; HHR23A; Rad23A
Species	Human
Gene ID	5886
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Rad23A (P54725).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54725
Protein Weight	40kDa
Shipping	Ice bag