


Anti-Rad23A Rabbit Polyclonal Antibody For Wb, Elisa - P54725

Item Number: CM0019618



Introduction

Anti-Rad23A rabbit polyclonal antibody validated for Western blotting (WB) and ELISA applications. Reacts with human and mouse samples; immunogen corresponds to full-length human Rad23A. Suitable for DNA repair and demethylation studies.

[Learn More](#)

Parameter	Value
Product Name	Rad23A Rabbit pAb
Remarks/Alias	HR23A; HHR23A; Rad23A
Species	Human
GeneID (Human)	5886
GeneID	5886
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-363 of human Rad23A (NP_005044.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54725
Protein Weight	40kDa
Shipping	Ice bag