


Anti-Rad23B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P54727

Item Number: CM0016451



Introduction

Rad23B rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Recognizes the 43kDa Rad23B protein, a key factor in DNA repair and histone demethylation.

[Learn More](#)

Parameter	Value
Product Name	Rad23B Rabbit pAb
Remarks/Alias	P58; HR23B; HHR23B; Rad23B
Species	Human
Gene ID	5887
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54727
Protein Weight	43kDa
Shipping	Ice bag