

Anti-Rad50 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92878

Item Number: CM0015985



Introduction

Anti-Rad50 polyclonal antibody for WB, IHC-P, ELISA; detects human and mouse RAD50, a key DNA repair protein. Suitable for cancer research and genomic stability studies.

[Learn More](#)

Parameter	Value
Product Name	Rad50 Rabbit pAb
Remarks/Alias	NBSLD; RAD502; hRad50; Rad50
Species	Human
GeneID (Human)	10111
GeneID	10111
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-270 of human Rad50 (NP_005723.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92878
Protein Weight	154kDa
Shipping	Ice bag