


Anti-Rad50 Rabbit Monoclonal Antibody For Wb, Elisa - Q92878

Item Number: CM0008155



Introduction

High-specificity rabbit monoclonal antibody targeting human RAD50, validated for WB and ELISA. Cross-reacts with mouse and rat RAD50. Ideal for DNA repair and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rad50 Rabbit mAb
Remarks/Alias	NBSLD; RAD502; hRad50; Rad50
Species	Human
Gene ID	10111
Immunogen	Synthetic peptide (aa 1200-1300 of human Rad50, Q92878)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92878
Protein Weight	154kDa
Shipping	Ice bag