


Anti-Rad50 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q92878

Item Number: CM0026386



Introduction

High-quality rabbit polyclonal antibody against human RAD50, validated for Western blot, IHC-P, IP, and ELISA. Reacts with human, mouse, and rat. Targets a key DNA repair protein in the MRN complex, essential for genomic stability and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rad50 Rabbit pAb
Remarks/Alias	NBSLD; RAD502; hRad50; Rad50
Species	Human
GeneID (Human)	10111
GeneID	10111
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Rad50 (NP_005723.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92878
Protein Weight	154kDa
Shipping	Ice bag