

Anti-Rad51 Rabbit Monoclonal Antibody For Wb, Elisa, Chip - Q06609

Item Number: CM0009998



Introduction

This rabbit monoclonal antibody targets DNA repair protein RAD51 homolog 1. Validated for WB, ELISA, and ChIP applications, it reacts with human, mouse, and rat samples, making it suitable for oncology and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	Rad51 Rabbit mAb
Remarks/Alias	RECA; BRCC5; FANCR; MRMV2; HRAD51; RAD51A; HsRad51; HsT16930; Rad51
Species	Human
GeneID (Human)	5888
GeneID	5888
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Rad51 (Q06609).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06609
Protein Weight	37kDa
Shipping	Ice bag