

Anti-Rad51 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q06609

Item Number: CM0008452



Introduction

High-quality Rabbit Monoclonal Antibody targeting human RAD51, validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Ideal for DNA repair, homologous recombination, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RAD51 Rabbit mAb
Remarks/Alias	RECA; BRCC5; FANCR; MRMV2; HRAD51; RAD51A; HsRad51; HsT16930
Species	Human
Gene ID	5888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human Rad51 (NP_002866.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06609
Protein Weight	37kDa
Shipping	Ice bag