

Anti-Rad51 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q06609

Item Number: CM0028453



Introduction

Rabbit polyclonal antibody against human DNA repair protein RAD51 (UniProt Q06609). Validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying homologous recombination and DNA repair mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rad51 Rabbit pAb
Remarks/Alias	RECA; BRCC5; FANCR; MRMV2; HRAD51; RAD51A; HsRad51; HsT16930; Rad51
Species	Human
Gene ID (Human)	5888
Gene ID	5888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human Rad51 (NP_002866.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06609
Protein Weight	37kDa
Shipping	Ice bag