

Anti-Rad51B Polyclonal Antibody For Wb And Elisa - O15315

Item Number: CM0025556



Introduction

RAD51B Rabbit Polyclonal Antibody validated for WB and ELISA, reacts with human, mouse, and rat. RAD51B is a DNA repair protein involved in homologous recombination, making this antibody ideal for DNA damage and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RAD51B Rabbit pAb
Remarks/Alias	REC2; R51H2; RAD51L1; RAD51B
Species	Human
GeneID (Human)	5890
GeneID	5890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-384 of human RAD51B (NP_598193.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15315
Protein Weight	42kDa
Shipping	Ice bag