

Anti-Rad51C Polyclonal Antibody For If/Icc And Elisa - O43502

Item Number: CM0016317



Introduction

RAD51C Rabbit pAb is a polyclonal antibody for IF/ICC and ELISA, targeting human RAD51C. Sourced from rabbit, it is ideal for investigating DNA repair, homologous recombination, and related pathways.

[Learn More](#)

Parameter	Value
Product Name	RAD51C Rabbit pAb
Remarks/Alias	FANCO; R51H3; BROVCA3; RAD51L2; RAD51C
Species	Human
Gene ID	5889
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RAD51C (NP_002867.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43502
Protein Weight	42kDa
Shipping	Ice bag