

Anti-Rad51C Rabbit Polyclonal Antibody For Wb, Elisa - O43502

Item Number: CM0029197



Introduction

Rabbit polyclonal antibody targeting human RAD51C, validated for Western blotting (WB) and ELISA. Detects ~42 kDa protein, cross-reacts with human and mouse. Ideal for DNA repair and recombination research.

[Learn More](#)

Parameter	Value
Product Name	RAD51C Rabbit pAb
Remarks/Alias	FANCO; R51H3; BROVCA3; RAD51L2; RAD51C
Species	Human
Gene ID (Human)	5889
Gene ID	5889
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human RAD51C (NP_002867.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43502
Protein Weight	42kDa
Shipping	Ice bag