

Anti-Rad51D Rabbit Polyclonal Antibody For Wb, Elisa - 075771

Item Number: CM0030176



Introduction

Anti-Rad51D rabbit polyclonal antibody for Western blot and ELISA applications. Validated in human samples, targeting the DNA repair protein RAD51 homolog 4 involved in homologous recombination repair.

[Learn More](#)

Parameter	Value
Product Name	Rad51D Rabbit pAb
Remarks/Alias	TRAD; R51H3; BROVCA4; RAD51L3; Rad51D
Species	Human
Gene ID (Human)	5892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-328 of human Rad51D (NP_002869.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075771
Protein Weight	35kDa
Shipping	Ice bag