


Anti-Rad52 Polyclonal Antibody For Wb, Elisa - P43351

Item Number: CM0027761



Introduction

Rabbit polyclonal antibody against human RAD52, suitable for Western blot and ELISA. Detects human, mouse, and rat RAD52. Ideal for DNA repair and recombination studies.

[Learn More](#)

Parameter	Specification
Product Name	Rad52 Rabbit pAb
Remarks/Alias	RAD52; Rad52
Species	Human
Gene ID (Human)	5893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-418 of human Rad52 (NP_602296.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43351
Protein Weight	46kDa
Shipping	Ice bag