

Anti-Rad54B Polyclonal Antibody For Wb, Ip, Elisa - Q9Y620

Item Number: CM0022852

Introduction

Rabbit polyclonal antibody targeting human RAD54B (103 kDa, SWISS: Q9Y620) for Western blot, immunoprecipitation, and ELISA applications. Cross-reacts with mouse. Ideal for DNA repair and homologous recombination research.

[Learn More](#)

Parameter	Value
Product Name	RAD54B Rabbit pAb
Remarks/Alias	RDH54; RAD54B
Species	Human
Gene ID (Human)	25788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 781-910 of human RAD54B (NP_036547.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y620
Protein Weight	103kDa
Shipping	Ice bag