

Anti-Rad54L2 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q9Y4B4

Item Number: CM0017143

Introduction

Rabbit polyclonal antibody against human RAD54L2 (ARIP4), validated for IHC-P, IF/ICC, ELISA. Reacts with human and rat RAD54L2, a helicase involved in androgen receptor transcription. Suitable for cancer research.

[Learn More](#)

Parameter	Value
Product Name	RAD54L2 Rabbit pAb
Remarks/Alias	ARIP4; HSPC325; SRISNF2L; RAD54L2
Species	Human
Gene ID (Human)	23132
Gene ID	23132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1178-1467 of human RAD54L2 (NP_055921.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9Y4B4
Protein Weight	163kDa
Shipping	Ice bag