


Anti-Rad9A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q99638

Item Number: CM0021218



Introduction

Rabbit polyclonal antibody against human RAD9A, validated in WB, IHC-P, and ELISA. Reacts with human and rat. Suitable for DNA repair and cell cycle checkpoint research. SWISS: Q99638.

[Learn More](#)

Parameter	Value
Product Name	RAD9A Rabbit pAb
Remarks/Alias	RAD9; RAD9A
Species	Human
GeneID (human)	5883
GeneID	5883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 162-391 of human RAD9A (NP_004575.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q99638
Protein Weight	43kDa
Shipping	Ice bag