

Anti-Radixin Polyclonal Antibody For Wb And Elisa - P35241

Item Number: CM0018173



Introduction

High-quality anti-Radixin rabbit polyclonal antibody for WB and ELISA applications. Detects human, mouse, and rat Radixin. Ideal for investigating actin cytoskeleton dynamics and membrane-cytoskeleton linkage.

[Learn More](#)

Parameter	Value
Product Name	Radixin Rabbit pAb
Remarks/Alias	DFNB24; Radixin
Species	Human
GeneID (Human)	5962
GeneID	5962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Radixin (NP_002897.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35241
Protein Weight	69kDa
Shipping	Ice bag