

Anti-Radixin Rabbit Monoclonal Antibody For Wb, Elisa - P35241

Item Number: CM0002421



Introduction

Anti-Radixin Rabbit Monoclonal Antibody targeting 69kDa protein, validated for WB and ELISA; cross-reacts with human, mouse, rat; for actin-membrane linkage and deafness studies.

[Learn More](#)

Parameter	Value
Product Name	Radixin Rabbit mAb
Remarks/Alias	DFNB24; Radixin
Species	Human
Gene ID (Human)	5962
Gene ID	5962
Immunogen	Synthetic peptide (amino acids 350-450 of human Radixin)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35241
Protein Weight	69kDa
Shipping	Ice bag