


Anti-Rae1 Rabbit Polyclonal Antibody For Wb, Elisa - P78406

Item Number: CM0030966



Introduction

Anti-RAE1 Rabbit Polyclonal Antibody targeting human RAE1 (mRNA export factor). Validated in WB and ELISA; cross-reacts with mouse and rat. Recombinant protein immunogen (1-368 aa). 41kDa. Ideal for nucleocytoplasmic transport and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	RAE1 Rabbit pAb
Remarks/Alias	Gle2, MIG14, MRNP41, Mnrp41, dj481F12.3, dj800J21.1, RAE1
Species	Human
Gene ID (Human)	8480
Gene ID	8480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-368 of human RAE1 (NP_003601.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78406
Protein Weight	41kDa
Shipping	Ice bag