

Anti-Raf1 Rabbit Monoclonal Antibody For Wb, Elisa - P04049

Item Number: CM0003525



Introduction

High-quality Anti-Raf1 Rabbit Monoclonal Antibody validated for WB and ELISA, with cross-reactivity to Human, Mouse, and Rat. Ideal for studying Raf1-mediated MAPK signaling pathways and cell fate regulation.

[Learn More](#)

Parameter	Value
Product Name	Raf1 Rabbit mAb
Remarks/Alias	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Raf1
Species	Human
GeneID	5894
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Raf1 (P04049)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04049
Protein Weight	73kDa
Shipping	Ice bag