

Anti-Raf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04049

Item Number: CM0017058



Introduction

Rabbit polyclonal antibody against human Raf1 (c-Raf). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Raf1 is a key kinase in the MAPK/ERK pathway, critical for cell proliferation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Raf1 Rabbit pAb
Remarks/Alias	NS5; CRAF; Raf-1; c-Raf; CMD1NN; Raf1
Species	Human
GeneID (Human)	5894
GeneID	5894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-350 of human Raf1 (NP_002871.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04049
Protein Weight	73kDa
Shipping	Ice bag