

Anti-Rag1 Rabbit Polyclonal Antibody For Wb, Elisa - P15918

Item Number: CM0013951

Introduction

This rabbit polyclonal antibody targets human RAG1 (V(D)J recombination-activating protein 1, also known as RNF74). Validated for Western blotting and ELISA, it cross-reacts with mouse and detects the ~119 kDa protein. Ideal for studies of lymphocyte development, immunodeficiency disorders, and V(D)J recombination mechanisms.

[Learn More](#)

Parameter	Value
Product Name	RAG1 Rabbit pAb
Remarks/Alias	RAG-1; RNF74; RAG1
Species	Human
Gene ID (Human)	5896
Gene ID	5896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human RAG1 (NP_000439.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P15918
Protein Weight	119kDa
Shipping	Ice bag