


Anti-Rag2 Rabbit Polyclonal Antibody For Wb, Elisa - P55895

Item Number: CM0016671



Introduction

Rabbit polyclonal antibody targeting human RAG2, a key V(D)J recombination protein. Validated for WB and ELISA with cross-reactivity to mouse. Ideal for lymphocyte development studies.

[Learn More](#)

Parameter	Value
Product Name	RAG2 Rabbit pAb
Remarks/Alias	RAG-2; RAG2
Species	Human
GeneID	5897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 258-527 of human RAG2 (NP_001230715.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P55895
Protein Weight	59kDa
Shipping	Ice bag