

Anti-Rai14 Rabbit Polyclonal Antibody For Wb And Elisa - Q9P0K7

Item Number: CM0017053

Introduction

Rabbit polyclonal antibody targeting human ankycorbin (RAI14). Suitable for Western blotting and ELISA with cross-reactivity to mouse and rat. Ideal for reproductive biology and actin regulation studies.

[Learn More](#)

Parameter	Value
Product Name	RAI14 Rabbit pAb
Remarks/Alias	RAI13, NORPEG, RAI14
Species	Human
GeneID	26064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-600 of human RAI14 (NP_056392.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9P0K7
Protein Weight	110kDa
Shipping	Ice bag