

Anti-Raidd/Cradd Rabbit Monoclonal Antibody For Wb, Elisa - P78560

Item Number: CM0010905



Introduction

High-specificity rabbit monoclonal antibody against human RAIDD/CRADD, validated in WB and ELISA. Suitable for apoptosis and TNFR-1 signaling pathway research. Recognizes human CRADD (23 kDa).

[Learn More](#)

Parameter	Value
Product Name	RAIDD/CRADD Rabbit mAb
Remarks/Alias	MRT34; RAIDD; RAIDD/CRADD
Species	Human
GeneID (Human)	8738
GeneID	8738
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-199 of human RAIDD/CRADD (NP_003796.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78560
Protein Weight	23 kDa
Shipping	Ice bag