

Anti-Raidd/Cradd Rabbit Polyclonal Antibody For Wb, Elisa - P78560

Item Number: CM0015603



Introduction

Human RAIDD/CRADD rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse RAIDD, a 23 kDa apoptosis adapter protein (SWISS P78560). Ideal for apoptosis and cell death studies.

[Learn More](#)

Parameter	Details
Product Name	RAIDD/CRADD Rabbit pAb
Remarks/Alias	MRT34; RAIDD; RAIDD/CRADD
Species	Human
GeneID (Human)	8738
GeneID	8738
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human RAIDD/CRADD (NP_003796.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78560
Protein Weight	23kDa
Shipping	Ice bag