

Anti-Rala Rabbit Polyclonal Antibody For Wb And Elisa - P11233

Item Number: CM0017130



Introduction

Anti-RALA rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat RALA (24 kDa). Suitable for research on GTPase signaling, cell migration, and oncogenic transformation.

[Learn More](#)

Parameter	Value
Product Name	RALA Rabbit pAb
Remarks/Alias	RAL; HINCONS; RALA
Species	Human
GeneID (Human)	5898
GeneID	5898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RALA (NP_005393.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11233
Protein Weight	24kDa
Shipping	Ice bag