

Anti-Ralbp1 Monoclonal Antibody For Ihc-P, If/Icc, Elisa - Q15311

Item Number: CM0012432



Introduction

High-quality rabbit monoclonal antibody against human RALBP1 (RLIP76), validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for Ral signaling and transport studies.

[Learn More](#)

Parameter	Value
Product Name	RALBP1 Rabbit mAb
Remarks/Alias	RIP1; RLIP1; RLIP76; RALBP1
Species	Human
GeneID (human)	10928
GeneID	10928
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human RALBP1 (Q15311).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15311
Protein Weight	76kDa
Shipping	Ice bag