

Anti-Ralbp1 Rabbit Polyclonal Antibody (Ko Validated) For Wb And Elisa - Q15311

Item Number: CM0015083



Introduction

Rabbit polyclonal antibody targeting RalA-binding protein 1 (RALBP1), also known as RLIP76. Validated for Western blot and ELISA, with cross-reactivity to human and mouse. Ideal for investigating Ral signaling, endocytosis, mitochondrial fission, and drug resistance mechanisms.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RALBP1 Rabbit pAb
Remarks/Alias	RIP1, RLIP1, RLIP76, P1
Species	Human
Gene ID (Human)	10928
Gene ID	10928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 396-655 of human RALBP1 (NP_006779.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15311
Protein Weight	76kDa
Shipping	Ice bag