

Anti-Ralgds Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q12967

Item Number: CM0021833



Introduction

Rabbit polyclonal antibody against human RALGDS (RalGEF), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying Ras/Ral signaling and bacterial clearance.

[Learn More](#)

Parameter	Value
Product Name	RALGDS Rabbit pAb
Remarks/Alias	RGF, RGDS, RalGEF, RALGDS
Species	Human
Gene ID (Human)	5900
Gene ID	5900
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 771-902 of human RALGDS (NP_001258705)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12967
Protein Weight	101kDa
Shipping	Ice bag