

Anti-Ralgps1 Rabbit Polyclonal Antibody For Wb And Elisa - Q5Js13

Item Number: CM0013669



Introduction

Rabbit polyclonal antibody targeting human RALGPS1 (RALGEF2) for western blot and ELISA. Suitable for studying RALA GTPase signaling, cytoskeletal dynamics, and membrane trafficking. Detects 62 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	RALGPS1 Rabbit pAb
Remarks/Alias	RALGEF2; RALGPS1A; RALGPS1
Species	Human
Gene ID (Human)	9649
Gene ID	9649
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RALGPS1 (NP_055451.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5J513
Protein Weight	62kDa
Shipping	Ice bag