


Anti-Raly Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukm9

Item Number: CM0013628



Introduction

RALY Rabbit pAb detects endogenous RALY protein in human, mouse, and rat samples. Validated for Western blot and ELISA. RALY is an RNA-binding protein involved in cholesterol biosynthesis and mRNA processing. UniProt Q9UKM9.

[Learn More](#)

Parameter	Value
Product Name	RALY Rabbit pAb
Remarks/Alias	P542; HNRPCL2; RALY
Species	Human
GeneID (human)	22913
GeneID	22913
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-306 of human RALY (NP_057951.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKM9
Protein Weight	32kDa
Shipping	Ice bag