

Anti-Ramp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60895

Item Number: CM0029250



Introduction

Anti-RAMP2 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Detects human, mouse, and rat RAMP2. Suitable for receptor complex studies and tissue localization.

[Learn More](#)

Parameter	Value
Product Name	RAMP2 Rabbit pAb
Remarks/Alias	RAMP2
Species	Human
GeneID (Human)	10266
GeneID	10266
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human RAMP2 (NP_005845.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60895
Protein Weight	20kDa
Shipping	Ice bag