

Anti-Ramp3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60896

Item Number: CM0029999



Introduction

RAMP3 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat samples. Investigates receptor activity-modifying protein 3 involvement in adrenomedullin and amylin receptor signaling.

[Learn More](#)

Parameter	Value
Product Name	RAMP3 Rabbit pAb
Remarks/Alias	RAMP3
Species	Human
Gene ID (Human)	10268
Gene ID	10268
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human RAMP3 (NP_005847.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60896
Protein Weight	17kDa
Shipping	Ice bag