


Anti-Ran Rabbit Monoclonal Antibody For Wb, Elisa - P62826

Item Number: CM0009324



Introduction

Rabbit monoclonal anti-Ran antibody for WB and ELISA applications. Detects Ran (24 kDa) in human, mouse, and rat. Ideal for studying nucleocytoplasmic transport and mitosis.

[Learn More](#)

Parameter	Value
Product Name	Ran Rabbit mAb
Remarks/Alias	TC4; Gsp1; ARA24; Ran
Species	Human
GeneID (Human)	5901
GeneID	5901
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 117-216 of human Ran (P62826).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62826
Protein Weight	24kDa
Shipping	Ice bag