


Anti-Ran Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62826

Item Number: CM0014931



Introduction

Rabbit polyclonal antibody targeting human Ran (GTP-binding nuclear protein), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying nucleocytoplasmic transport and mitosis.

[Learn More](#)

Parameter	Value
Product Name	Ran Rabbit pAb
Remarks/Alias	TC4; Gsp1; ARA24; Ran
Species	Human
GeneID (Human)	5901
GeneID	5901
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-216 of human Ran (NP_006316.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62826
Protein Weight	24 kDa
Shipping	Ice bag